

# JD GUIDED SURGERY WORKFLOW

JDentalCare accompanies you all along the way and puts at your disposal the experience and capabilities of JDLab, a laboratory of the latest generation to which you can turn

1

## THE DENTIST SENDS TO THE JD LAB:

- CT scan in DICOM format
- Diagnostic wax-up copy
- Templates
- Medical prescription
- Evobite

2

## THE JD LAB PREPARES:

- Model scans with Evobite
- Guided-surgery digital planning
- Creation of the surgical guide

3

## THE JD LAB SENDS TO DENTIST:

- CT scan in DICOM format
- Diagnostic wax-up copy
- Templates
- Medical Prescription
- Evobite
- +
- Surgical guide
- Chewing index
- 3D printed model
- Digital planning summary

